

Date TK

Contact: Liz Dowling, (760) 942-2544
Dowling & Dennis Public Relations
E-mail: liz@dowlingdennis.net

NEW PRODUCT RELEASE

SwabCap from Excelsior Medical Disinfects Luer Access Valve Cap

Excelsior Medical has launched the SwabCap™ Luer Access Valve Disinfection Cap, designed to disinfect swabable luer access valves prior to line access. SwabCap twists directly onto the valve, passively disinfecting the valve top and threads without activating the luer access valve. The simple application and removal process promotes technique standardization and caregiver compliance with luer access valve cleaning protocols, including the new (January 2010) Joint Commission requirement.

SwabCap also protects patients by acting as a physical barrier to touch and airborne contamination, before line access and between line accesses. The cap provides up to 96 hours of protection under normal conditions if not removed.

-- END --