

The **only** saline flush syringe that offers a **SwabCap®** built into the syringe plunger.

SwabFlush™



SwabFlush™

Offers clinicians the convenience of having SwabCap **there when they need it**, after the final flush!

When to use SwabFlush


- When a saline flush/lock is required to finish a patient's IV therapy
- Use SwabFlush to administer the final saline:
 - After catheter insertion
 - After medication delivery
 - After blood is withdrawn or delivered through the catheter

SwabCap®

**NON-VENTING
DISINFECTION CAP
CONTAINS 70% IPA**

Single Use Only
Sterile Packaging
Luer Lock Design

Simple to use!



PRODUCT DESCRIPTION

SwabFlush™ 10/10mL Saline Flush Syringe
SwabFlush™ 5/10mL Saline Flush Syringe
SwabFlush™ 3/10mL Saline Flush Syringe

Code

E0103-01
E0103-55
E0103-35

Units/Box

100
100
100

Units/Case

600
600
800



www.excelsiormedical.com

Call toll free 800-487-4276 or visit www.excelsiormedical.com to find out more about SwabCap®

SwabCap®

Valve disinfection is a critical part of IV catheter maintenance and infection prevention.

Let SwabCap® help you:

- Achieve disinfection and protection of needleless connectors
- Standardize disinfection technique and drive compliance to protocol
- Visually monitor and measure compliance
- Save nursing time

**NON-VENTING
DISINFECTION CAP
CONTAINS 70% IPA**

Single Use Only
Sterile Packaging
Luer Lock Design



SwabCap. . . Designed to meet your Highest Quality Standards

Sterile Packaging

No cross contamination fears for your patient

Thread Cover Seal

Alcohol stays in and does not dry out

Medical Grade Sponge Pad


No particulates upon applying or removal

Soft and Low Profile

Gentle on patient's skin

Always at your fingertips

Integrated into your flush syringe



Patent pending **thread cover design** helps ensure complete bathing of the connector's top and threads

SwabFlush™



The **best** of both technologies



Corporate Headquarters
Excelsior Medical Corporation
1933 Heck Avenue, Neptune, NJ 07753

Toll Free 800-487-4276
Main Phone 732-776-7525
Fax 732-776-7600

Web: www.excelsiormedical.com
E-mail: info@excelsiormedical.com